



Rest assured with NiteOwl, your companion for anytime anywhere sleep therapy.

Airline travel • Battery backup • Camping

SPECIFICATIONS

Cell Type: Rechargeable Lithium-ion

Capacity: 99.9 Wh

Output Voltage: 12V-24V DC (depending on DC output cable)

DC Input Voltage Range: 10V-28V

Quick-charge 3.0 USB Output Voltage: 5V | 2.4A (maximum)

Charge Input: 10V-28V

Charge Time: 2-3 hrs @ 90W

*charge times are per battery and using the ResMed Air 10 AC power supply

Warranty: 1 year



Ask about NiteOwl at your local sleep therapy equipment provider or contact us.
We can help you find a retailer in your area.

www.niteowlcpapbattery.com • 866-718-0557 • sleepniteowl@gmail.com



Rest assured with NiteOwl, your companion for anytime anywhere sleep therapy.

Use as a Battery Backup for Uninterrupted Therapy During Power Outages



Pictured: Single battery backup setup with ResMed Air 11 Series

Power Your PAP Devices On the Go



Pictured: Dual battery standard setup with ResMed Air 11 Series

Ask about NiteOwl at your local sleep therapy equipment provider or contact us.
We can help you find a retailer in your area.

www.niteowlcpapbattery.com • 866-718-0557 • sleepniteowl@gmail.com



Rest assured with NiteOwl, your companion for anytime anywhere sleep therapy.



NITEOWL-1-KIT

Universal battery kit. For use with PR System One 60s, PR Dreamstation 1&2, RESMED A10, A11, and Airmini.



NITEOWL-2-KIT

Universal Dual (2) battery kit. For use with PR System One 60s, PR Dreamstation 1 & 2, RESMED A10, A11, and Airmini.



NITEOWL-BATT

Spare battery only, with Bridge Cable.



PHRESP-CAB

Cable only for use with 12Vdc Philips Respironics System One 60s, and Dreamstation 1 and 2.

Ask about NiteOwl at your local sleep therapy equipment provider or contact us.
We can help you find a retailer in your area.

www.niteowlcpapbattery.com • 866-718-0557 • sleepniteowl@gmail.com

RESMED-CABKIT

Cable kit only for use with 24Vdc ResMed Airsense 10 & 11 series, and Airmini travel device. Includes Power Supply Adapters. *included in the RESMED-CABKIT



A10-CAB*

Airsense 10 output cable only



AM-A11-CAB*

Airmini and Airsense 11 output cable only



A11-PSA*

Airsense 11 Power Supply Adapter



AM-PSA*

Airmini Power Supply Adapter

Ask about NiteOwl at your local sleep therapy equipment provider or contact us.
We can help you find a retailer in your area.



Rest assured with NiteOwl, your companion for anytime anywhere sleep therapy.



BRIDGE-CAB

Bridge cable only. Used for connecting to NiteOwl batteries in dual mode.



AUTO-CAB

Mobile Power Adapter "car charger" for charging via 12VDC outlet.



AC-CHARGER

AC Charger. 100-240VAC wall charging

Ask about NiteOwl at your local sleep therapy equipment provider or contact us.
We can help you find a retailer in your area.

www.niteowlcpapbattery.com • 866-718-0557 • sleepniteowl@gmail.com